



EDITOR'S CHOICE

This economical diode laser features plenty of power and a touchscreen interface ideal for new laser users.

**A** Diode laser



**B** Disposable prophy tip



**C** Interproximal reduction system



**D** Restoration alignment system



**E** Air drying system

**AMD Lasers LLC**

**A Powerful, economical diode laser**

Featuring a touch screen interface and up to 7W of power, the **Picasso** diode laser is designed as an economical option for dentists looking to add laser services to an operator. With continuous and pulse modes, eight customizable presets and an adjustable aiming beam, the laser is suited to a variety of procedures. The unit features multiple menu language options, is compatible with U.S. and international power outlets and ships with one fiber, one handpiece, three pairs of goggles, a training DVD, a "Quick Start" manual and a 3-year warranty.

866-999-2635 • [amdlasers.com](http://amdlasers.com)

**Groman Inc.**

**B Disposable air polish prophy**

Designed to simultaneously provide full mouth removal of extrinsic stains, dental plaque and soft debris while polishing, the **ProphyMaster™** is a waterless, disposable air polishing prophylaxis tip integrated with a high-speed evacuator tube. The tube is said to capture a significant portion of the aerated powder dust as it ricochets off the tooth surface. Each tip combines the suctioning power of the chair-side high-speed evacuation and low-speed air supply to provide minimal overspray in a single-handed operation; flow control is provided via the chair foot.

800-906-9298  
[gromandental.com](http://gromandental.com)

**Spee Corp.**

**C Fast, safe interproximal reduction**

Created to help prevent reducing too much enamel during an interproximal reduction, the **Spee Slenderizer System** also addresses clinician concerns about patient movement when the diamond is in close proximity to soft tissue. The system's 16:1 Straight Reduction Handpiece slows down the rotation of the diamond disk without losing any of the power needed to properly perform IPR procedures.

866-369-1978 • [speecorp.com](http://speecorp.com)

**Common Sense Dental**

**D Create centered restorations**

Indicated for use on veneers, anterior implant placements, dentures and to verify bonding, **PerfectAlign** is designed to help improve dentists' communications with labs to create centered restorations. It is said to allow the dentist to clearly see the entirety of anterior teeth and communicate details such as the patient's midline, vertical axial inclinations, the horizontal plane and to verify the incisal edge position of the central incisors.

888-853-5773  
[commonsensedental.com](http://commonsensedental.com)

**Air Techniques Inc.**

**E Non-stop membrane dryer system**

Able to precipitate humidity through a membrane, the **membrane dryer system** presents an efficient solution for dental air drying. The compressed air is cooled to ambient temperature by the aluminum ribs of a cooling element. Air then flows through 1,200 thin membrane fibers that dehumidify it. On the other side of the membrane fibers, dry scavenging air absorbs the humidity again—said to precipitate three times more water than conventional systems.

800-AIR-TECH (800-247-8324)  
[airtechniques.com](http://airtechniques.com)